



Homemade Microbial Media

By Kenneth Rainis, Adjunct Professor, Genesee Community College

MATERIALS

Containers: tin-foil muffin tins, clear deli containers (containers must be heat-resistant)

Measuring cup (heat resistant)

Mixing bowl / container (microwavable)

Aluminum foil / plastic wrap

Spoon

Gelatin (Knox brand or equivalent)

Water

Bouillon cube (beef)

Table sugar

PROCEDURE

1. Mix together 250ml (1 cup) of water with 2 packets (~12g) gelatin into a mixing container.
2. Mix in 1 bouillon cube.
3. Mix in 10g sugar (~10ml or 2tsp)
4. Stir to dissolve/mix all contents thoroughly
5. Bring the mixture to a boil in a microwave by repeated heating (use 1-minute intervals), and watching carefully until the gelatin (or agar) is completely dissolved.
6. Let the mixture cool for 15 minutes. (WARNING: **HOT LIQUID!**)
7. Carefully pour the hot medium into clean containers until 1/2 or 1/3 full.
8. Carefully cover (wrap) each container with either aluminum foil or plastic wrap.
9. Allow to cool and harden. Do not agitate the containers while cooling/hardening is taking place.
10. Use the microbe food **within 3 days**. Store in a refrigerator or other cool, dry location.

NOTES

1) Gelatin will melt if it gets too warm (depending upon the grade) ~ 35 ° C). Some microbe strains will liquefy gelatin (that's why agar agar is used in formulations). Agar Agar – can be found in most Asian grocery stores or baking supply stores & is the better choice of a solidification medium.

2) One envelope of powdered gelatin, ~7g is ~35ml or ~2.5 tsp

3) 5g (~1 tsp) bouillon granules can be used in place of cube